



Dr. perfect V-line


PDO THREAD

Features

Only Dr. Perfect V-Line

Has Cavern

Having ZERO pitch tunnel shape
Strongly resist press by integumentary system
Produce collagen not only outside but also inside of thread
Spin-less insert



Has Multi Thread

Insert 14 threads at one treatment.

Caps opener & Station

You can open five caps of needle at one try with mounted caps opener. And put opened products into the holes of station, and then operator can pick them up one by one.



CE Certificate

From the first issue date (Jan. 10th .2014)
CE committee has strengthened its audit process very critically recently. All the CE certificate holders must pass post assessment every year.
Many of manufactures have lost their CE certificate and trying to get new certificate from other notified body.
So just few of manufactures maintain its original CE certificate.
It proves that how well they control quality of products.

Blister Packing (Made of Tyvek)

Tyvek can protect the contents from humidity and air much more than any other material.
Blister packing enables automatic production.

Drying process

Packing and sealing in vacuum chamber with 0 % humidity.

Expiry Date

Due to Tyvek packing and drying process, we guarantee 2 years expiry date.



CAVERN

Patent no. 10-2016-0010527

NEEDLE GAUGE	NEEDLE LENGTH	USP	THREAD LENGTH	THREAD LENGTH (STRETCHED)	USE
29G	50mm	6-0	20mm	271mm	Eye wrinkles
27G	50mm	5-0	20mm	254mm	Eye bag
25G	60mm	5-0	20mm	289mm	Glabella
23G(CN)	60mm	4-0	20mm	273mm	Sunken Temple filling
23G(CN)	60mm	4-0	30mm	398mm	Sunken Cheek filling Sunken Temple filling Nasolabial Jawline
21G(CN)	60mm	3-0	25mm	370mm	Midcheek groove

- Strongly resist press by integumentary System due to very tight tunnel shape
- Powerful volume-up
- Much more collagen collect inside and outside
- Spin-less insert
- Filling sunken cheeks

MULTI THREAD

Patent no. 10-1487230

NEEDLE GAUGE	NEEDLE LENGTH	USP	THREAD LENGTH	USE
23G(CN)	38mm	5-0 × 14	25mm	Nasolabial
	60mm	5-0 × 14	45mm	

- Insert 14 threads at one treatment



Forehead wrinkles

Glabella

Crow's feet


Tear trough

Midcheek groove

Nasolabial


Marionette line

Chin Crease





COG(C)



Needle type / Cannula type

NEEDLE GAUGE	NEEDLE LENGTH	THREAD LENGTH	USP
23G	60 / 90	105 / 210	(3-0)
21G	60 / 90	140 / 170	(2-0)
19G	100	150 / 210	(0-0)

CLEOKO(CLK)

Nasal treatment



NEEDLE GAUGE	NEEDLE LENGTH	THREAD LENGTH	USP
21G	60mm	90mm	(0-0)
19G	38mm	28mm	(0-0)
19G(Cog)	38mm	46mm	(0-0)

MONO(N)




NEEDLE GAUGE	NEEDLE LENGTH	THREAD LENGTH	USP
30G	25 / 38	30 / 50	6-0
29G	25 / 38 / 50 / 60	30 / 50 / 70 / 90	6-0
27G	38 / 50 / 60	50 / 70 / 90	5-0
26G	38 / 50 / 60 / 90	50 / 70 / 90 / 150	5-0

DOUBLE(D)




NEEDLE GAUGE	NEEDLE LENGTH	THREAD LENGTH	USP
29G	25 / 38 / 50 / 60	30 / 50 / 70 / 90	(7-0)×2
27G	38 / 50 / 60	50 / 70 / 90	(6-0)×2
26G	38 / 50 / 60 / 90	50 / 70 / 90 / 150	(7-0)+(6-0)
25G	60 / 90	90 / 150	(6-0)×2

SCREW(S)



NEEDLE GAUGE	NEEDLE LENGTH	THREAD LENGTH	USP
30G	25 / 38	30 / 50	6-0
29G	25 / 38 / 50 / 60	30 / 50 / 70 / 90	6-0
27G	38 / 50 / 60	50 / 70 / 90	5-0
26G	38 / 50 / 60 / 90	50 / 70 / 90 / 150	5-0

DOUBLE SCREW(DS)



NEEDLE GAUGE	NEEDLE LENGTH	THREAD LENGTH	USP
29G	25 / 38 / 50 / 60	30 / 50 / 70 / 90	(7-0)×2
27G	38 / 50 / 60	50 / 70 / 90	(7-0)×2
26G	38 / 50 / 60 / 90	50 / 70 / 90 / 150	(7-0)+(6-0)
25G	60 / 90	90 / 150	(6-0)×2



Dr. perfect V-line

PACKAGE



CAVERN, COG,
MULTI THREAD
100EA / Individual box



MONO, DOUBLE,
SCREW, DOUBLE SCREW
200EA / Individual box



CARTON
2,400EA / CTN

CE 2195 ISO 13485 GMP KFDA

neoDr NEO Dr. INC.
www.neodr.net

Tel : +82-70-8766-2569, Fax : +82-33-734-1863

E-mail : sales@neodr.com, <http://www.neodr.net>

#2-204, 205, 207, Medical Industry Complex, 42-10 Taejanggongdan-gil, Wonju-city,
Gangwon-do, Republic of Korea 26311