



Reyoungel®

Cross-linked Hyaluronic Acid- Dermal Filler



See What
We Can Do
For Your
Beauty!



Reyoungel Cross-linked Hyaluronic Acid- Dermal Filler

What is hyaluronic Acid ?

Hyaluronic acid is a polysaccharide and is one of very few substances that it is identical in all species and tissue types. Hyaluronic acid also exists naturally in our bodies and its most important function is to carry and bind water. It plays an important role giving volume to the skin, shape to the eyes and elasticity to the joints. It can be found in the eyes, tendons, articulations, and above all in the dermis of the skin. With the years our skin produces less hyaluronic acid, which is the result of normal aging process. Your skin becomes thinner and as a result, it becomes loose and starts to wrinkle.

What is Reyoungel ?

Reyoungel is biodegradable gel made of non-animal cross-linked hyaluronic acid. They offer a complete solution for treating facial lines, wrinkles, and lips. It can also be used in facial sculpting by enhancing natural facial contours and features. There are three types of products can be Chosen: Reyoungel Derm, Reyoungel Area and Reyoungel Lips. These three products make it possible to quickly and effectively meet the vast majority of patient's needs for correction and augmentation of wrinkles and folds.



Reyoungel Products info

	Reyoungel Derm	Reyoungel Area	Reyoungel Lips	Reyoungel Derm	Reyoungel Lips
Concentration	20mg/ml stabilized Hyaluronic acid	20mg/ml stabilized Hyaluronic acid	20mg/ml stabilized Hyaluronic acid	20mg/ml stabilized Hyaluronic acid	20mg/ml stabilized Hyaluronic acid
Approximate number of gel particles/ml	10000	10000	10000	10000	10000
Pharmaceutical grade (mg/ml)	0.30-0.40	0.15-0.30	0.30-0.40	0.3-1.34	0.30-0.33
Recommended Indications	Deep facial wrinkles and folds, like horizontal forehead, glabellar, nasolabial, marionette, and perioral wrinkles, and smile lines	Medium wrinkles, such as nasolabial and nasolabial folds, nasolabial and perioral wrinkles, and smile lines	Thin superficial lines, such as nasolabial, perioral, and periorbital wrinkles, and smile lines	Lip line definition and lip volume enhancement, and smile lines	The best definition and lip volume enhancement, and smile lines
Method of Injekt	Deep layer of injection and/or subcutaneous of injection	Medium layer of injection	Deep layer of injection	Deep layer of injection and/or subcutaneous of injection	Deep layer of injection and/or subcutaneous of injection
Volume of syringe	1ml/ml	1ml/ml	1ml/ml	1ml/ml	1ml/ml
Needle size	30G	30G	30G	30G	30G
Duration	6-12months	6-12months	6-12months	6-12months	6-12months

Reyoungel Benefits of Reyoungel

Safe & Natural :

- * From non-animal origin and high-purity
- * Excellent biodegradability in the skin
- * No side effects, high quality, less invasive
- * Fill the wrinkles and smooth your skin without inhibiting facial expression.

Instant result :

- * The result lasts for 18months. A touch up treatment after a year needs the half of the product.
- * Reyoungel contains the high concentration of stabilized Hyaluronic acid, which remains to longer residence time.

Fast & Effective :

- * Non-surgical within 30 minutes procedure
- * No hospital stay needed, short recovery
- * Noticeable results after the first treatment.
- * Subtle breast re-shaping and fullness



Reyoungel Technical Data

Active Component	Sodium Hyaluronate (Cross-linked)
Concentration	20 mg/ml
Volume	1.0ml/3ml
pH	6.8-7.2
Package	one Pre-filled syringe with two 30G needles
Storage	Between 2-8°C protected from sunlight and frost



Reyoungel Certifications

Free Sale Certificate
European Co. Mark
Et Design - Examination Certificate
Et ISO 13485: 2012

