



Ultra Pure Water Purification System Model Pure-12

[Get Quote](#)[Download Catalog](#)

SPECIFICATION

FEATURES

Output of Pure Water	12 L/hr
Output of Ultra Pure Water	>1.5 L/min
Resistivity of ultra pure water (25°C)	18.2 MΩ-cm
TOC	10 ppb
Bacteria	<0.1 CFU/ml
Particle (>0.1µm)	<1/ml
Heavy Metal Ion	<0.1 ppb
Conductivity of RO water	1-5µs/cm*
Dimension	550x390x610 mm
Power Supply	110-220V 50/60 Hz
Weight	20 Kg
Catalog No.	461101248